

Sample Planchets

Models L-7525-371 and L-7525-371-01



Ludlum Measurements, Inc.



Part Numbers:
Aluminum L-7525-371 and
Stainless Steel L-7525-371-01

Introduction

Ludlum Measurements, Inc manufactures sample planchets for use with its sample counters, or those of other manufacturers that accommodate a standard 5.1 cm (2 inch) planchet. The devices are available in either aluminum or stainless steel in quantities of not less than 500 pieces.

Specifications

INDICATED USE:

SIZE: 0.37 x 5.08 x 5.08 cm (0.15 x 2.0 x 2.0 in.) (H x W x L)

MATERIAL CHOICE:

Aluminum Model L-7525-371

Stainless Steel Model L-7525-371-01