Model 3-98 with 44-3 Alpha Beta Gamma Survey Meter



Features

- Dual Purpose Detection
- . 125 Alpha Beta & Gamma Contamination
- 4-Range Analog Ratemeter
- 0-500 kcpm
- CE Certified
- Rugged
- Complete Turn-Key System

Specifications

Part Number: 48-1135

INDICATED USE: alpha, beta, low-energy gamma survey

DETECTORS: controlled by selector switch

INTERNAL DETECTOR:

Type: pancake-style GM detector with $1.7 \pm 0.3 \text{ mg/cm}^2 \text{ window,}$ stainless steel protective screen (79% open)

Area:

Active: 15 cm² Open: 12 cm²

Efficiency $(4\pi \text{ geometry}): 5\%^{-14}\text{C}; 23\%^{-90}\text{Sr}/^{90}\text{Y}; 19\%^{-99}\text{Tc}; 33\%^{-32}\text{P}; 15\%^{-239}\text{Pu}$

Sensitivity: typically 3300 cpm/mR/hr (137Cs gamma)
Alpha-beta Shield: retractable 211 mg/cm² phenolic slide

EXTERNAL DETECTOR: Model 44-3

Type: low-energy gamma scintillator

Efficiency (4π geometry): 33.5%–¹²⁵I (based on ¹²⁹I efficiency of 18%)

Background (10 μ R/hr): < 350 cpm Sensitivity: 675 cpm/ μ R/hr (125|)

Energy Range: recommended 5 keV–60 keV **CONNECTOR**: series "C" (others available)

METER DIAL: 0-5 kcpm, BAT TEST (others available)

MULTIPLIERS: x0.1, x1, x10, x100

LINEARITY: reading within ±10% of true value with detector connected

AUDIO: built-in unimorph speaker with ON/OFF switch (greater than 60 dB at 0.61 m [2 ft]) **CALIBRATION CONTROLS**: accessible from front of instrument (protective cover provided)

HIGH VOLTAGE: adjustable from 400–1500 volts

THRESHOLD: -30 mV ±10 mV

RESPONSE: toggle switch for FAST (4 seconds) or SLOW (22 seconds) from 10% to 90% of final reading

RESET: pushbutton to zero meter

POWER: 2 "D" cell batteries (housed in sealed compartment that is externally accessible)

BATTERY LIFE: typically 600 hours with alkaline batteries, battery condition can be checked on meter

METER: 6.4 cm (2.5 in.) arc, 1 mA analog type

CONSTRUCTION: cast and drawn aluminum with beige powder coat

TEMPERATURE RANGE: -20 to 50 °C (-4 to 122 °F)

SIZE: 22.9 x 8.9 x 21.6 cm (9 x 3.5 x 8.5 in.) (H x W x L), with handle

WEIGHT: 2.5 kg (5.6 lb) including batteries and detector





View of bottom of instrument showing retractable slide.

