

# Model L-301

## Table Top Densitometer



Ludlum Measurements, Inc.

### Introduction

The Model L-301 Table Top and Model L-331C Portable Densitometers are easy to use precision instruments designed to provide highly accurate and repeatable (black & white) optical density readings. This makes them an ideal tool for Processor Quality Assurance.

The readings provided by the Model L-301 (and L-331C) will alert you to fluctuations in processing conditions and allow you to take the necessary corrective action before film quality becomes an issue. The Model L-301 offers an optional RS-232 interface.

### Specifications


Measuring range: 0-5.0 D with 2 and 3 mm apertures; 0-4.0 D with 1 mm aperture



Part Number: 99-9600

# Model L-331C

## Portable Densitometer



Part Number: 99-9601

The Model L-331C Portable Densitometer is an easy-to-use precision instrument designed to provide highly accurate and repeatable (black and white) optical density readings. This makes it an ideal tool for processor quality assurance. The readings provided by the Model L-331C (and the Model L-301) will alert you to fluctuations in processing conditions and allow you to take the necessary corrective action before film quality becomes an issue.

**MEASURING RANGE:** 0 to 5.0 D with 2 and 3mm apertures; 0 to 4.0 D with 1 mm aperture  
A 200 Vdc version is available upon request for the same price.