

Model L-929-100 & Model L-929-092 Shielded Syringe Holder & Carrier

Model L-929-100

Lead Shielded Syringe Holder

Part Number: 99-9301

The lead shielded syringe holder will accommodate unshielded, isotope-containing syringes ranging from 1 to 20 cc. The holder has 1.3 cm (0.5 in.) thick lead shielding, and an extra-wide base prevents accidental tipping.

Specifications

OUTSIDE DIMENSIONS (H x W): 16.5 x 9.8 cm (6.5 x 3.9 in.)

INSIDE DIMENSIONS (H x D): 14.3 x 2 cm (5.6 x 0.8 in.)

SHIELDING: 1.3 cm (0.5 in.) lead

WEIGHT: 2.7 kg (6 lb)



Model L-929-092

Shielded Stainless Steel Syringe Carrier

Part Number: 99-9300

The Model L-929-092 Stainless Steel Syringe Carrier is constructed entirely of stainless steel and provides shielding equivalent to 3.175 mm (1/8 in.) lead. The syringe carrier will store and transport syringes, vials or ampules up to 22.9 cm (9 in.) long. The lid is overlapping with a latch to prevent streaming and opens 180° for easy access. There are four rubber feet attached to the bottom of the carrier.

Specifications

OUTSIDE DIMENSIONS WITH HANDLE (H x W x L):

12.3 x 9.9 x 24.8 cm (4.83 x 3.91 x 9.75 in.)

INSIDE DIMENSIONS (H x W x L):

5.4 x 5.7 x 22.9 cm (2.125 x 2.25 x 9 in.)

MATERIAL: 304 stainless steel

SHIELDING: 9.525 mm (3/8 in.) thick type 304 stainless steel (thicker material used on cover)

LEAD EQUIVALENCY: 3.175 mm (1/8 in.)

WEIGHT: about 4 kg (8 lb)

