

# Model L-DeCon-1 Decontamination Kit

*Radiation Detection for a Safer World*



Ludlum Measurements, Inc.



Part Number 4525-561

## Introduction

The LUDLUM Decontamination Kit contains the essential equipment, labeling materials, protective clothing, decontamination solvents and cleaning supplies, needed to deal with a basic radioactive spill or routine decontamination problem that may be encountered in the Nuclear Medicine Department or radiologic research laboratory. The 30 gallon drum serves both as a container for the kit components as well as a waste transfer and temporary storage (decay) vessel.

## Standard Kit Components

- 1 30 Gallon Drum ( hazmat rated )
- 2 each Disposable Coverall and Shoe Covers
- 2 each N-95 Mask
- 2 pairs Disposable Gloves
- 1 Radiacwash Spray 1 liter \*\*
- 1 Box Radiacwash Towlettes (100/box)
- 10 each 30 Gallon Trash Bags
- 1 each 18" Crucible Tongs
- 2 each Mop heads
- 1 each Mop handle
- 1 Plastic (3 Gal) Pail
- 1 scrub brush ( 6")
- 1 Roll ( 3 " x 1000' ) "Caution Rad Area " Barricade Tape
- 2 Each Laminated (7" x 10 " ) "Radioactive Material" sign
- 3 Each vinyl (5" x 7") " Radioactive Material" Labels
- 1 Roll Absorbent Paper Towels

## Options

- L-Decon-2 (Deluxe Kit - Part Number: 4525-562) which includes: DeCon Gel and additional Coveralls, Gloves and Shoe Covers
- DeCon Gel / Spray Bottle, 1 liter Bottle - Part Number 03-5648