

Model 180-15-3 Sample Holder

Part Number: 47-3865



Specifications

INDICATED USE: repeatable geometry counting of wipes, filter paper, or slides

SUGGESTED DETECTOR: Model 43-92, Model 43-93, Model 44-142

SAMPLE SIZE: 5.1 x 1.5 cm (2 x 0.6 in.) (Dia x L) maximum

SIZE: 2.5 x 9.4 x 18.5 cm (1 x 3.7 x 7.3 in.) (H x W x L)

WEIGHT: 1.1 kg (2.3 lb)