


Model 180-15 Sample Holder

Radiation Detection for a Safer World



Ludlum Measurements, Inc.



Part Number: 47-1111

Specifications

INDICATED USE: repeatable geometry counting of wipes, filter paper, or slides

SUGGESTED DETECTOR: Model 44-40

SAMPLE SIZE: 5.1 x 1.5 cm (2 x 0.6 in.) (Dia x L) maximum

SIZE: 2.5 x 9.4 x 18.5 cm (1 x 3.7 x 7.3 in.) (H x W x L)

WEIGHT: 1.1 kg (2.3 lb)

Sample Planchets

Choose aluminum or stainless steel version:

ALUMINUM: Model L-7525-371

STAINLESS STEEL: Model L-7525-371-01

DIAMETER: 5.1 cm (2.0 in.)

THICKNESS: 3.2 mm (0.125 in.)

Minimum order quantity: 500

