

# Model 180-15 Sample Holder



Ludlum Measurements, Inc.



## Specifications

Part Number: 47-1111

**INDICATED USE:** repeatable geometry counting of wipes, filter paper, or slides

**SUGGESTED DETECTOR:** Model 44-40

**SAMPLE SIZE:** 5.1 x 1.5 cm ( 2 x 0.6 in.) (Dia x L) maximum

**SIZE:** 2.5 x 9.4 x 18.5 cm (1 x 3.7 x 7.3 in.) (H x W x L)

**WEIGHT:** 1.1 kg (2.3 lb)

## Sample Planchets

Choose aluminum or stainless steel version:

ALUMINUM: Model L-7525-371

STAINLESS STEEL: Model L-7525-371-01

STANDARD SIZES (DIAMETER X THICKNESS):

5.1 x 1.5 cm (2.0 x 0.6 in.)

5.7 x 1.3 cm (2.25 x 0.5 in.)

5.1 x 0.32 cm (2.0 x 0.125 in.)

2.5 x 0.32 cm (1.0 x 0.125 in.)

Other sizes available upon request.



Minimum order quantity: 500

Ludlum Measurements, Inc. P.O. Box 810, Sweetwater, Texas 79556

**Web:** <http://www.ludlums.com> **Tel:** 800-622-0828 / 325-235-5494 / **Fax:** 325-235-4672 / **Email:** [sales@ludlums.com](mailto:sales@ludlums.com)

Note: specifications subject to change without notification. We are not responsible for errors or omissions.

June 2024