

# Model 180-28

## Sample Holder for Model 26-Series Friskers



Ludlum Measurements, Inc.



Left: Model 180-28  
Right: Model 180-28 with Model 26

### Specifications

Part Number: 47-3948

**INDICATED USE:** repeatable geometry counting of wipes, filter paper, or slides

**SAMPLE SIZE:** 5.1 cm (2 in.) diameter maximum

**SAMPLE TRAY:** selectable spacing of 0.32 cm (0.125 in.), 0.64 cm (0.25 in.), 1.3 cm (0.5 in.), 2.5 cm (1 in.), and 5.1 cm (2 in.) from detector

**SUGGESTED DETECTOR:** Model 26 Integrated Frisker

**CONSTRUCTION:** anodized aluminum frame and sample tray

**SIZE:** 9.1 x 7.6 x 15.2 cm (3.6 x 3 x 6 in.) (H x W x L)

**WEIGHT:** 0.40 kg (0.9 lb)

### Model 26-Series Integrated Friskers

The Model 180-28 Sample Holder is designed to accommodate the design of the Model 26-Series Integrated Friskers to make testing for contamination on wipes, filter paper, or slides simple and effective. As shown in the two views at right, the frisker head rests in the specially designed top of the sample holder while the handle fits in the cutout. Model 26 is available separately.

**Model 26** Part Number: 48-3885

