

Model 43-78

Alpha Sample Counter



Ludlum Measurements, Inc.



Part Number: 47-2180

Specifications

INDICATED USE: large sample alpha counting
SUGGESTED INSTRUMENTS: Models 2000, 2200
SAMPLE SIZE (Maximum): 12.7 x 0.46 cm (5 x 0.18 in.) (Dia x L)
SAMPLE HOLDER: anodized aluminum tray
SHIELDING: none
SCINTILLATOR: ZnS(Ag) adhered to a 0.6 cm (0.25 in.) light pipe
PHOTOMULTIPLIER TUBE: 12.7 cm (5 in.) diameter
OPERATING VOLTAGE: 500–1200 volts
WINDOW: windowless
EFFICIENCY (4 π): 37%–²³⁰Th; 37%–²³⁹Pu
BACKGROUND: alpha: 3 cpm or less
SIZE: 19.1 x 31.8 cm (7.5 x 12.6 in.) (Dia x L)
WEIGHT: 6.8 kg (15 lb)

Available:

Model	Sample Size (Dia x Depth)	Part Number
Model 43-78	12.7 x 0.46 cm (5 x 0.46 in.)	47-2180
Model 43-78-1	12.7 x 0.46 cm (5 x 0.46 in.)	47-2307
Model 43-78-2	7.6 x 0.46 cm (3 x 0.18 in.)	47-2620
Model 43-78-3	12.7 x 0.46 cm (5 x 0.46 in.)	47-2762

Ludlum Measurements, Inc. P.O. Box 810, Sweetwater, Texas 79556

Web: <http://www.ludlums.com> Tel: 800-622-0828 / 325-235-5494 / Fax: 325-235-4672 / Email: ludlum@ludlums.com

Note: specifications subject to change without notification. We are not responsible for errors or omissions.

March 2015